




FABP6 Antibody

Product Code	CSB-PA007955ESR1HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P51161
Immunogen	Recombinant Human Gastrotropin protein (1-128AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Gastrotropin (GT) (Fatty acid-binding protein 6) (Ileal lipid-binding protein) (ILBP) (Intestinal 15 kDa protein) (I-15P) (Intestinal bile acid-binding protein) (I-BABP), FABP6, ILBP ILLBP
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	FABP6

Image



Immunohistochemistry of paraffin-embedded human colon cancer using CSB-PA007955ESR1HU at dilution of 1:100

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.