




NCBP2 Antibody

Product Code	CSB-PA015521ESR2HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P52298
Immunogen	Recombinant Human Nuclear cap-binding protein subunit 2 protein (1-156AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Nuclear cap-binding protein subunit 2 (20 kDa nuclear cap-binding protein) (Cell proliferation-inducing gene 55 protein) (NCBP 20 kDa subunit) (CBP20) (NCBP-interacting protein 1) (NIP1), NCBP2, CBP20
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	NCBP2

Image



Immunohistochemistry of paraffin-embedded human colon cancer using CSB-PA015521ESR2HU at dilution of 1:100

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.