




ACTR2 Antibody

Product Code	CSB-PA001248ESR1HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P61160
Immunogen	Recombinant Human Actin-related protein 2 protein (1-160AA)
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	ELISA, WB, IHC; Recommended dilution: WB:1:500-1:1000, IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Actin-related protein 2 (Actin-like protein 2), ACTR2, ARP2
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	ACTR2

Image



Western blot

All lanes: Actin-related protein 2 antibody at 3µg/ml

Lane 1: Mouse spleen tissue

Lane 2: Jurkat whole cell lysate


Lane 3: MCF-7 whole cell lysate

Secondary

Goat polyclonal to rabbit IgG at 1/10000 dilution

Predicted band size: 45, 46 kDa

Observed band size: 45 kDa



Immunohistochemistry of paraffin-embedded human tonsil tissue using CSB-PA001248ESR1HU at dilution of 1:100



Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.