




FADD Antibody

Product Code	CSB-PA623788ESR1HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q13158
Immunogen	Recombinant Human FAS-associated death domain protein (1-208AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	FAS-associated death domain protein (FAS-associating death domain-containing protein) (Growth-inhibiting gene 3 protein) (Mediator of receptor induced toxicity) (Protein FADD), FADD, MORT1
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	FADD

Image


Immunohistochemistry of paraffin-embedded human melanoma using CSB-PA623788ESR1HU at dilution of 1:100

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.