





MAP2K5 Antibody

Product Code	CSB-PA615667ESR1HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q13163
Immunogen	Recombinant Human Dual specificity mitogen-activated protein kinase kinase 5 protein (1-190AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Dual specificity mitogen-activated protein kinase kinase 5 (MAP kinase kinase 5) (MAPKK 5) (EC 2.7.12.2) (MAPK/ERK kinase 5) (MEK 5), MAP2K5, MEK5 MKK5 PRKMK5
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	MAP2K5

Image



Immunohistochemistry of paraffin-embedded human testis tissue using CSB-PA615667ESR1HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human gallbladder tissue using CSB-PA615667ESR1HU at dilution of 1:100



Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.