



BLMH Antibody

Product Code	CSB-PA002716ESR1HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q13867
Immunogen	Recombinant Human Bleomycin hydrolase protein (1-300AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Bleomycin hydrolase (BH) (BLM hydrolase) (BMH) (EC 3.4.22.40), BLMH
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	BLMH
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.