





ELOA Antibody

Product Code	CSB-PA619867ESR2HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q14241
Immunogen	Recombinant Human Elongin-A protein (679-798AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Elongin-A (EloA) (Elongin 110 kDa subunit) (RNA polymerase II transcription factor SIII subunit A1) (SIII p110) (Transcription elongation factor B polypeptide 3), ELOA, TCEB3
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	ELOA

Image



Immunohistochemistry of paraffin-embedded human tonsil tissue using CSB-PA619867ESR2HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human colon cancer using CSB-PA619867ESR2HU at dilution of 1:100



Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.