





# SMC1A Antibody

|                            |                                                                                                                                             |
|----------------------------|---------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Product Code</b>        | CSB-PA614804ESR1HU                                                                                                                          |
| <b>Storage</b>             | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.                                                                               |
| <b>Uniprot No.</b>         | Q14683                                                                                                                                      |
| <b>Immunogen</b>           | Recombinant Human Structural maintenance of chromosomes protein 1A protein (1-130AA)                                                        |
| <b>Raised In</b>           | Rabbit                                                                                                                                      |
| <b>Species Reactivity</b>  | Human                                                                                                                                       |
| <b>Tested Applications</b> | ELISA, IHC; Recommended dilution: IHC:1:20-1:200                                                                                            |
| <b>Form</b>                | Liquid                                                                                                                                      |
| <b>Conjugate</b>           | Non-conjugated                                                                                                                              |
| <b>Storage Buffer</b>      | PBS with 0.02% sodium azide, 50% glycerol, pH7.3.                                                                                           |
| <b>Purification Method</b> | Antigen Affinity Purified                                                                                                                   |
| <b>Isotype</b>             | IgG                                                                                                                                         |
| <b>Clonality</b>           | Polyclonal                                                                                                                                  |
| <b>Alias</b>               | Structural maintenance of chromosomes protein 1A (SMC protein 1A) (SMC-1-alpha) (SMC-1A) (Sb1.8), SMC1A, DXS423E KIAA0178 SB1.8 SMC1 SMC1L1 |
| <b>Immunogen Species</b>   | Homo sapiens (Human)                                                                                                                        |
| <b>Research Area</b>       | Cell Biology                                                                                                                                |
| <b>Target Names</b>        | SMC1A                                                                                                                                       |

## Image



Immunohistochemistry of paraffin-embedded human tonsil tissue using CSB-PA614804ESR1HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human colon cancer using CSB-PA614804ESR1HU at dilution of 1:100