




TDG Antibody

Product Code	CSB-PA623816ESR2HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q13569
Immunogen	Recombinant Human G/T mismatch-specific thymine DNA glycosylase protein (141-410AA)
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	ELISA, WB; Recommended dilution: WB:1:500-1:2000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity purified
Isotype	IgG
Clonality	Polyclonal
Alias	G/T mismatch-specific thymine DNA glycosylase (EC 3.2.2.29) (Thymine-DNA glycosylase) (hTDG), TDG
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	TDG

Image



Western blot

All lanes: G/T mismatch-specific thymine DNA glycosylase antibody at 8µg/ml + Mouse brain tissue

Secondary

Goat polyclonal to rabbit IgG at 1/10000 dilution

Predicted band size: 46 kDa

Observed band size: 46 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.