





# CST1 Antibody

<b>Product Code</b>	CSB-PA006088ESR1HU
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	P01037
<b>Immunogen</b>	Recombinant Human Cystatin-SN protein (20-141AA)
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Human
<b>Tested Applications</b>	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
<b>Form</b>	Liquid
<b>Conjugate</b>	Non-conjugated
<b>Storage Buffer</b>	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
<b>Purification Method</b>	Antigen Affinity Purified
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Alias</b>	Cystatin-SN (Cystatin-SA-I) (Cystatin-1) (Salivary cystatin-SA-1), CST1
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Research Area</b>	Cell Biology
<b>Target Names</b>	CST1

## Image



Immunohistochemistry of paraffin-embedded human prostate tissue using CSB-PA006088ESR1HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human salivary gland tissue using CSB-PA006088ESR1HU at dilution of 1:100