




RPS4X Antibody

Product Code	CSB-PA020456ESR2HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P62701
Immunogen	Recombinant Human 40S ribosomal protein S4, X isoform protein (66-263AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity purified
Isotype	IgG
Clonality	Polyclonal
Alias	40S ribosomal protein S4, X isoform (SCR10) (Single copy abundant mRNA protein) (Small ribosomal subunit protein eS4), RPS4X, CCG2 RPS4 SCAR
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	RPS4X

Image



Immunohistochemistry of paraffin-embedded human colon cancer using CSB-PA020456ESR2HU at dilution of 1:100

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.