





CALD1 Antibody

Product Code	CSB-PA004441DSR1HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q05682
Immunogen	Recombinant Human Caldesmon protein (1-240AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Caldesmon (CDM), CALD1, CAD CDM
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	CALD1


Image



Immunohistochemistry of paraffin-embedded human prostate tissue using CSB-PA004441DSR1HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human colon tissue using CSB-PA004441DSR1HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human colon cancer using CSB-PA004441DSR1HU at dilution of 1:100

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.