





MAPK7 Antibody

Product Code	CSB-PA013465ESR1HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q13164
Immunogen	Recombinant Human Mitogen-activated protein kinase 7 protein (1-310AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC, IF; Recommended dilution: IHC:1:20-1:200, IF:1:50-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Mitogen-activated protein kinase 7 (MAP kinase 7) (MAPK 7) (EC 2.7.11.24) (Big MAP kinase 1) (BMK-1) (Extracellular signal-regulated kinase 5) (ERK-5), MAPK7, BMK1 ERK5 PRKM7
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	MAPK7

Image



Immunohistochemistry of paraffin-embedded human kidney tissue using CSB-PA013465ESR1HU at dilution of 1:100



Immunofluorescent analysis of HeLa cells using CSB-PA013465ESR1HU at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)



Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.