





ETV6 Antibody

Product Code	CSB-PA007860ESR1HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P41212
Immunogen	Recombinant Human Transcription factor ETV6 protein (1-300AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity purified
Isotype	IgG
Clonality	Polyclonal
Alias	Transcription factor ETV6 (ETS translocation variant 6) (ETS-related protein Tel1) (Tel), ETV6, TEL TEL1
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	ETV6


Image



Immunohistochemistry of paraffin-embedded human breast cancer using CSB-PA007860ESR1HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human skin tissue using CSB-PA007860ESR1HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human epityphlon tissue using CSB-PA007860ESR1HU at dilution of 1:100