



ST14 Antibody

Product Code	CSB-PA022746ESR2HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9Y5Y6
Immunogen	Recombinant Human Suppressor of tumorigenicity 14 protein (566-855AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Suppressor of tumorigenicity 14 protein (EC 3.4.21.109) (Matriptase) (Membrane-type serine protease 1) (MT-SP1) (Prostamin) (Serine protease 14) (Serine protease TADG-15) (Tumor-associated differentially-expressed gene 15 protein), ST14, PRSS14 SNC19 TADG15
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	ST14
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.