




IRAK3 Antibody

Product Code	CSB-PA011811ESR1HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9Y616
Immunogen	Recombinant Human Interleukin-1 receptor-associated kinase 3 protein (367-596AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Interleukin-1 receptor-associated kinase 3 (IRAK-3) (EC 2.7.11.1) (IL-1 receptor-associated kinase M) (IRAK-M), IRAK3
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	IRAK3

Image



Immunohistochemistry of paraffin-embedded human epityphlon tissue using CSB-PA011811ESR1HU at dilution of 1:100

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.