




ANTXR1 Antibody

Product Code	CSB-PA001832ESR1HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9H6X2
Immunogen	Recombinant Human Anthrax toxin receptor 1 protein (315-564AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity purified
Isotype	IgG
Clonality	Polyclonal
Alias	Anthrax toxin receptor 1 (Tumor endothelial marker 8), ANTXR1, ATR TEM8
Immunogen Species	Homo sapiens (Human)
Research Area	Tags & Cell Markers
Target Names	ANTXR1

Image


Immunohistochemistry of paraffin-embedded human lung cancer using CSB-PA001832ESR1HU at dilution of 1:100

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.