




BMF Antibody

| | |
|----------------------------|---|
| Product Code | CSB-PA850383ESR1HU |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | Q96LC9 |
| Immunogen | Recombinant Human Bcl-2-modifying factor protein (1-184AA) |
| Raised In | Rabbit |
| Species Reactivity | Human |
| Tested Applications | ELISA, IHC; Recommended dilution: IHC:1:20-1:200 |
| Form | Liquid |
| Conjugate | Non-conjugated |
| Storage Buffer | PBS with 0.02% sodium azide, 50% glycerol, pH7.3. |
| Purification Method | Antigen Affinity Purified |
| Isotype | IgG |
| Clonality | Polyclonal |
| Alias | Bcl-2-modifying factor, BMF |
| Immunogen Species | Homo sapiens (Human) |
| Research Area | Cell Biology |
| Target Names | BMF |

Image



Immunohistochemistry of paraffin-embedded human liver cancer using CSB-PA850383ESR1HU at dilution of 1:100