



GDF2 Antibody

Product Code	CSB-PA892325GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9UK05
Immunogen	Human GDF2
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,IHC
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.1% sodium azide and 50% glycerol pH 7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	GDF-2,Bone morphogenetic protein 9,BMP-9,GDF2,BMP9
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	GDF2