



# GDF2 Antibody

<b>Product Code</b>	CSB-PA892325GA01HU
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	Q9UK05
<b>Immunogen</b>	Human GDF2
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Human,Mouse,Rat
<b>Tested Applications</b>	ELISA,IHC
<b>Form</b>	Liquid
<b>Conjugate</b>	Non-conjugated
<b>Storage Buffer</b>	PBS with 0.1% sodium azide and 50% glycerol pH 7.3.
<b>Purification Method</b>	Antigen Affinity Purified
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Alias</b>	GDF-2,Bone morphogenetic protein 9,BMP-9,GDF2,BMP9
<b>Product Type</b>	Polyclonal Antibody
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Target Names</b>	GDF2
<b>Usage</b>	For Research Use Only. Not for use in diagnostic or therapeutic procedures.