



PDX1 Antibody

Product Code	CSB-PA017746GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P52945
Immunogen	Human PDX1
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB,IHC
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.1% sodium azide and 50% glycerol pH 7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	PDX-1,Glucose-sensitive factor,GSF,Insulin promoter factor 1,IPF-1,Insulin upstream factor 1,IUF-1,Islet/duodenum homeobox-1, IDX-1,Somatostatin-transactivating factor 1,STF-1,PDX1,IPF1, STF1
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	PDX1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.