




TREX1 Antibody

Product Code	CSB-PA865133ESR2HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9NSU2
Immunogen	Recombinant Human Three-prime repair exonuclease 1 protein (1-270AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Three-prime repair exonuclease 1 (EC 3.1.11.2) (3'-5' exonuclease TREX1) (Deoxyribonuclease III) (DNase III), TREX1
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	TREX1

Image



Immunohistochemistry of paraffin-embedded human small intestine tissue using CSB-PA865133ESR2HU at dilution of 1:100