





PADI4 Antibody

Product Code	CSB-PA890757ESR2HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9UM07
Immunogen	Recombinant Human Protein-arginine deiminase type-4 protein (1-250AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Protein-arginine deiminase type-4 (EC 3.5.3.15) (HL-60 PAD) (Peptidylarginine deiminase IV) (Protein-arginine deiminase type IV), PADI4, PAD4 PADI5 PDI5
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	PADI4

Image



Immunohistochemistry of paraffin-embedded human breast cancer using CSB-PA890757ESR2HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human spleen tissue using CSB-PA890757ESR2HU at dilution of 1:100

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.