





TOMM20 Antibody

Product Code	CSB-PA618983ESR2HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q15388
Immunogen	Recombinant Human Mitochondrial import receptor subunit TOM20 homolog protein (20-145AA)
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	ELISA, WB, IHC, IF; Recommended dilution: WB:1:500-1:2000, IHC:1:20-1:200, IF:1:50-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Mitochondrial import receptor subunit TOM20 homolog (Mitochondrial 20 kDa outer membrane protein) (Outer mitochondrial membrane receptor Tom20), TOMM20, KIAA0016
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	TOMM20


Image



Western blot
 All lanes: Mitochondrial import receptor subunit TOM20 homolog antibody at 3μg/ml
 Lane 1: Mouse brain tissue
 Lane 2: Hela whole cell lysate
 Lane 3: Mouse liver tissue
 Secondary
 Goat polyclonal to rabbit IgG at 1/10000 dilution
 Predicted band size: 16 kDa
 Observed band size: 16 kDa



Immunohistochemistry of paraffin-embedded human kidney tissue using CSB-PA618983ESR2HU at dilution of 1:100



Immunofluorescent analysis of HepG2 cells using CSB-PA618983ESR2HU at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)