





CD48 Antibody

| | |
|----------------------------|---|
| Product Code | CSB-PA004941ESR1HU |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | P09326 |
| Immunogen | Recombinant Human CD48 antigen protein (27-125AA) |
| Raised In | Rabbit |
| Species Reactivity | Human |
| Tested Applications | ELISA, IHC; Recommended dilution: IHC:1:20-1:200 |
| Form | Liquid |
| Conjugate | Non-conjugated |
| Storage Buffer | PBS with 0.02% sodium azide, 50% glycerol, pH7.3. |
| Purification Method | Antigen Affinity Purified |
| Isotype | IgG |
| Clonality | Polyclonal |
| Alias | CD48 antigen (B-lymphocyte activation marker BLAST-1) (BCM1 surface antigen) (Leukocyte antigen MEM-102) (SLAM family member 2) (SLAMF2) (Signaling lymphocytic activation molecule 2) (TCT.1) (CD antigen CD48), CD48, BCM1 BLAST1 |
| Immunogen Species | Homo sapiens (Human) |
| Research Area | Immunology |
| Target Names | CD48 |

Image



Immunohistochemistry of paraffin-embedded human colon cancer using CSB-PA004941ESR1HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human epityphlon tissue using CSB-PA004941ESR1HU at dilution of 1:100