



BIRC8 Antibody

Product Code	CSB-PA846680ESR2HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96P09
Immunogen	Recombinant Human Baculoviral IAP repeat-containing protein 8 protein (1-236AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Baculoviral IAP repeat-containing protein 8 (Inhibitor of apoptosis-like protein 2) (IAP-like protein 2) (ILP-2) (Testis-specific inhibitor of apoptosis), BIRC8, ILP2
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	BIRC8
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.