




# RAB5A Antibody

<b>Product Code</b>	CSB-PA019213ESR1HU
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	P20339
<b>Immunogen</b>	Recombinant Human Ras-related protein Rab-5A protein (1-215AA)
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Human, Mouse
<b>Tested Applications</b>	ELISA, WB, IHC; Recommended dilution: WB:1:500-1:5000, IHC:1:20-1:200
<b>Form</b>	Liquid
<b>Conjugate</b>	Non-conjugated
<b>Storage Buffer</b>	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
<b>Purification Method</b>	Antigen Affinity Purified
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Alias</b>	Ras-related protein Rab-5A, RAB5A, RAB5
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Research Area</b>	Neuroscience
<b>Target Names</b>	RAB5A

## Image



### Western blot

All lanes: Ras-related protein Rab-5A antibody at 5µg/ml

Lane 1: Mouse brain tissue

Lane 2: HL60 whole cell lysate


Lane 3: Hela whole cell lysate

Secondary

Goat polyclonal to rabbit IgG at 1/10000 dilution

Predicted band size: 24, 23 kDa

Observed band size: 24 kDa



Immunohistochemistry of paraffin-embedded human colon cancer using CSB-PA019213ESR1HU at dilution of 1:100



---

**Usage**

For Research Use Only. Not for use in diagnostic or therapeutic procedures.