





RAB5A Antibody

Product Code	CSB-PA019213ESR1HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P20339
Immunogen	Recombinant Human Ras-related protein Rab-5A protein (1-215AA)
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	ELISA, WB, IHC; Recommended dilution: WB:1:500-1:5000, IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Ras-related protein Rab-5A, RAB5A, RAB5
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	RAB5A

Image



Western blot
 All lanes: Ras-related protein Rab-5A antibody at 5µg/ml
 Lane 1: Mouse brain tissue
 Lane 2: HL60 whole cell lysate
 Lane 3: Hela whole cell lysate
 Secondary
 Goat polyclonal to rabbit IgG at 1/10000 dilution
 Predicted band size: 24, 23 kDa
 Observed band size: 24 kDa



Immunohistochemistry of paraffin-embedded human colon cancer using CSB-PA019213ESR1HU at dilution of 1:100