





ATG16L1 Antibody

Product Code	CSB-PA734569ESR1HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q676U5
Immunogen	Recombinant Human Autophagy-related protein 16-1 protein (90-300AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB, IHC; Recommended dilution: WB:1:200-1:2000, IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Autophagy-related protein 16-1 (APG16-like 1), ATG16L1, APG16L
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	ATG16L1


Image



Western blot
 All lanes: ATG16L1 antibody at 2µg/ml + MCF-7 whole cell lysate
 Secondary
 Goat polyclonal to rabbit IgG at 1/10000 dilution
 Predicted band size: 69, 67, 54, 48, 50 kDa
 Observed band size: 69 kDa



Immunohistochemistry of paraffin-embedded human placenta tissue using CSB-PA734569ESR1HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human spleen tissue using CSB-PA734569ESR1HU at dilution of 1:100

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.