





CSRNP1 Antibody

Product Code	CSB-PA822303ESR2HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96S65
Immunogen	Recombinant Human Cysteine/serine-rich nuclear protein 1 protein (1-200AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Cysteine/serine-rich nuclear protein 1 (CSRNP-1) (Axin-1 up-regulated gene 1 protein) (Protein URAX1) (TGF-beta-induced apoptosis protein 3) (TAIP-3), CSRNP1, AXUD1 TAIP3
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	CSRNP1

Image



Immunohistochemistry of paraffin-embedded human colon cancer using CSB-PA822303ESR2HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human placenta tissue using CSB-PA822303ESR2HU at dilution of 1:100



Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.