




PLIN2 Antibody

Product Code	CSB-PA860767ESR1HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q99541
Immunogen	Recombinant Human Perilipin-2 protein (168-437AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB, IHC; Recommended dilution: WB:1:200-1:1000, IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Perilipin-2 (Adipophilin) (Adipose differentiation-related protein) (ADRP), PLIN2, ADFP
Immunogen Species	Homo sapiens (Human)
Research Area	Cardiovascular
Target Names	PLIN2

Image



Western blot

All lanes: PLIN2 antibody at 4µg/ml

Lane 1: HepG2 whole cell lysate


Lane 2: Hela whole cell lysate

Secondary

Goat polyclonal to rabbit IgG at 1/10000 dilution

Predicted band size: 48 kDa

Observed band size: 48 kDa



Immunohistochemistry of paraffin-embedded

human breast cancer using CSB-

PA860767ESR1HU at dilution of 1:100



Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.