




GRK2 Antibody

Product Code	CSB-PA001394DSR1HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P25098
Immunogen	Recombinant Human Beta-adrenergic receptor kinase 1 protein (440-689AA)
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	ELISA, WB; Recommended dilution: WB:1:200-1:1000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Beta-adrenergic receptor kinase 1 (Beta-ARK-1) (EC 2.7.11.15) (G-protein coupled receptor kinase 2), GRK2, ADRBK1 BARK BARK1
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	GRK2

Image



Western blot

All lanes: ADRBK1 antibody at 3µg/ml + Mouse brain tissue

Secondary

Goat polyclonal to rabbit IgG at 1/10000 dilution

Predicted band size: 80 kDa

Observed band size: 80 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.