



ALKBH3 Antibody

Product Code	CSB-PA847209ESR1HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96Q83
Immunogen	Recombinant Human Alpha-ketoglutarate-dependent dioxygenase alkB homolog 3 protein (1-286AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Alpha-ketoglutarate-dependent dioxygenase alkB homolog 3 (EC 1.14.11.54) (Alkylated DNA repair protein alkB homolog 3) (hABH3) (DEPC-1) (Prostate cancer antigen 1), ALKBH3, ABH3 DEPC1
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	ALKBH3
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.