




# ALKBH8 Antibody

<b>Product Code</b>	CSB-PA842644ESR1HU
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	Q96BT7
<b>Immunogen</b>	Recombinant Human Alkylated DNA repair protein alkB homolog 8 protein (1-230AA)
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Human
<b>Tested Applications</b>	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
<b>Form</b>	Liquid
<b>Conjugate</b>	Non-conjugated
<b>Storage Buffer</b>	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
<b>Purification Method</b>	Antigen Affinity Purified
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Alias</b>	Alkylated DNA repair protein alkB homolog 8 (EC 1.14.11.-) (Probable alpha-ketoglutarate-dependent dioxygenase ABH8) (S-adenosyl-L-methionine-dependent tRNA methyltransferase ABH8) (tRNA (carboxymethyluridine(34)-5-O)-methyltransferase ABH8) (EC 2.1.1.229), ALKBH8, ABH8
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Research Area</b>	Epigenetics and Nuclear Signaling
<b>Target Names</b>	ALKBH8

## Image



Immunohistochemistry of paraffin-embedded human colon cancer using CSB-PA842644ESR1HU at dilution of 1:100

## Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.