



ALKBH8 Antibody

Product Code	CSB-PA842644ESR2HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96BT7
Immunogen	Recombinant Human Alkylated DNA repair protein alkB homolog 8 protein (1-230AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity purified
Isotype	IgG
Clonality	Polyclonal
Alias	Alkylated DNA repair protein alkB homolog 8 (EC 1.14.11.-) (Probable alpha-ketoglutarate-dependent dioxygenase ABH8) (S-adenosyl-L-methionine-dependent tRNA methyltransferase ABH8) (tRNA (carboxymethyluridine(34)-5-O)-methyltransferase ABH8) (EC 2.1.1.229), ALKBH8, ABH8
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	ALKBH8
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.