




AMY1A Antibody

Product Code	CSB-PA356351ESR1HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P04745
Immunogen	Recombinant Human Alpha-amylase 1 protein (16-200AA)
Raised In	Rabbit
Species Reactivity	Human, Rat
Tested Applications	ELISA, WB; Recommended dilution: WB:1:500-1:2000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Alpha-amylase 1 (EC 3.2.1.1) (1,4-alpha-D-glucan glucanohydrolase 1) (Salivary alpha-amylase), AMY1A; AMY1B; AMY1C, AMY1; AMY1; AMY1
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	AMY1A

Image



Western blot
 All lanes: AMY1A antibody at 7µg/ml + Rat heart tissue
 Secondary
 Goat polyclonal to rabbit IgG at 1/10000 dilution
 Predicted band size: 58 kDa
 Observed band size: 58 kDa

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.