




AMY1A Antibody

Product Code	CSB-PA356351ESR2HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P04745
Immunogen	Recombinant Human Alpha-amylase 1 protein (16-200AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Alpha-amylase 1 (EC 3.2.1.1) (1,4-alpha-D-glucan glucanohydrolase 1) (Salivary alpha-amylase), AMY1A; AMY1B; AMY1C, AMY1; AMY1; AMY1
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	AMY1A

Image



Immunohistochemistry of paraffin-embedded human salivary gland tissue using CSB-PA356351ESR2HU at dilution of 1:100