



Sialic acid-binding lectin Antibody

Product Code	CSB-PA321644LA01RJL
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P18839
Immunogen	Recombinant Rana japonica Sialic acid-binding lectin protein (1-111AA)
Raised In	Rabbit
Species Reactivity	Rana japonica
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Sialic acid-binding lectin (EC 3.1.27.-)
Immunogen Species	Rana japonica (Japanese reddish frog)
Research Area	Others
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.