



Sialic acid-binding lectin Antibody, HRP conjugated

| | |
|---------------------|--|
| Product Code | CSB-PA321644LB01RJL |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | P18839 |
| Immunogen | Recombinant Rana japonica Sialic acid-binding lectin protein (1-111AA) |
| Raised In | Rabbit |
| Species Reactivity | Rana japonica |
| Tested Applications | ELISA |
| Form | Liquid |
| Conjugate | HRP |
| Storage Buffer | Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4 |
| Purification Method | >95%, Protein G purified |
| Isotype | IgG |
| Clonality | Polyclonal |
| Alias | Sialic acid-binding lectin (EC 3.1.27.-) |
| Immunogen Species | Rana japonica (Japanese reddish frog) |
| Research Area | Others |