




CASP14 Antibody

Product Code	CSB-PA004546ESR1HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P31944
Immunogen	Recombinant Human Caspase-14 protein (1-242AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Caspase-14 (CASP-14) (EC 3.4.22.-) [Cleaved into: Caspase-14 subunit p17, mature form; Caspase-14 subunit p10, mature form; Caspase-14 subunit p20, intermediate form; Caspase-14 subunit p8, intermediate form], CASP14
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	CASP14

Image



Immunohistochemistry of paraffin-embedded human skin tissue using CSB-PA004546ESR1HU at dilution of 1:100

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.