





CDK9 Antibody

Product Code	CSB-PA005079ESR1HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P50750
Immunogen	Recombinant Human Cyclin-dependent kinase 9 protein (153-372AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Cyclin-dependent kinase 9 (EC 2.7.11.22) (EC 2.7.11.23) (C-2K) (Cell division cycle 2-like protein kinase 4) (Cell division protein kinase 9) (Serine/threonine-protein kinase PITALRE) (Tat-associated kinase complex catalytic subunit), CDK9, CDC2L4 TAK
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	CDK9

Image



Immunohistochemistry of paraffin-embedded human endometrial cancer using CSB-PA005079ESR1HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human colon cancer using CSB-PA005079ESR1HU at dilution of 1:100