





CEBPG Antibody

Product Code	CSB-PA005184DSR1HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P53567
Immunogen	Recombinant Human CCAAT/enhancer-binding protein gamma protein (1-150AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	CCAAT/enhancer-binding protein gamma (C/EBP gamma), CEBPG
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	CEBPG

Image



Immunohistochemistry of paraffin-embedded human tonsil tissue using CSB-PA005184DSR1HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human placenta tissue using CSB-PA005184DSR1HU at dilution of 1:100