





CNOT8 Antibody

Product Code	CSB-PA871398ESR2HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9UFF9
Immunogen	Recombinant Human CCR4-NOT transcription complex subunit 8 protein (1-292AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	CCR4-NOT transcription complex subunit 8 (EC 3.1.13.4) (CAF1-like protein) (CALIFp) (CAF2) (CCR4-associated factor 8) (Caf1b), CNOT8, CALIF POP2
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	CNOT8

Image



Immunohistochemistry of paraffin-embedded human prostate cancer using CSB-PA871398ESR2HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human colon cancer using CSB-PA871398ESR2HU at dilution of 1:100



Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.