



IolA Antibody, FITC conjugated

| | |
|----------------------------|--|
| Product Code | CSB-PA419575LC01EJF |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | A7ZYJ5 |
| Immunogen | Recombinant Escherichia coli O9:H4 Outer-membrane lipoprotein carrier protein (22-203AA) |
| Raised In | Rabbit |
| Species Reactivity | Escherichia coli O9:H4 |
| Form | Liquid |
| Conjugate | FITC |
| Storage Buffer | Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4 |
| Purification Method | >95%, Protein G purified |
| Isotype | IgG |
| Clonality | Polyclonal |
| Alias | Outer-membrane lipoprotein carrier protein, IolA |
| Immunogen Species | Escherichia coli O9:H4 (strain HS) |
| Research Area | Others |
| Target Names | IolA |
| Usage | For Research Use Only. Not for use in diagnostic or therapeutic procedures. |