



lolA Antibody, Biotin conjugated

Product Code	CSB-PA419575LD01EJF
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	A7ZYJ5
Immunogen	Recombinant Escherichia coli O9:H4 Outer-membrane lipoprotein carrier protein (22-203AA)
Raised In	Rabbit
Species Reactivity	Escherichia coli O9:H4
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Outer-membrane lipoprotein carrier protein, lolA
Immunogen Species	Escherichia coli O9:H4 (strain HS)
Research Area	Others
Target Names	lolA
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.