



L2 Antibody, HRP conjugated

Product Code	CSB-PA356939LB01HMN
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P06793
Immunogen	Recombinant Human papillomavirus type 18 Minor capsid protein L2 protein (1-462AA)
Raised In	Rabbit
Species Reactivity	Human papillomavirus type 18
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Minor capsid protein L2, L2
Immunogen Species	Human papillomavirus type 18
Research Area	Others
Target Names	L2