



L2 Antibody, HRP conjugated

Product Code	CSB-PA356939LB01HMN
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P06793
Immunogen	Recombinant Human papillomavirus type 18 Minor capsid protein L2 protein (1-462AA)
Raised In	Rabbit
Species Reactivity	Human papillomavirus type 18
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Minor capsid protein L2, L2
Immunogen Species	Human papillomavirus type 18
Research Area	Others
Target Names	L2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.