



TGFB1 Antibody, Biotin conjugated

Product Code	CSB-PA023446LD01BO
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P18341
Immunogen	Recombinant Bovine Transforming growth factor beta-1 proprotein protein (279-390AA)
Raised In	Rabbit
Species Reactivity	Bovine
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Transforming growth factor beta-1 proprotein [Cleaved into: Latency-associated peptide (LAP)], (Transforming growth factor beta-1) (TGF-beta-1) TGFB1
Immunogen Species	Bos taurus (Bovine)
Research Area	Others
Target Names	TGFB1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.