



# ampC Antibody, HRP conjugated

|                     |                                                                                                                   |
|---------------------|-------------------------------------------------------------------------------------------------------------------|
| Product Code        | CSB-PA326492HB01EZX                                                                                               |
| Storage             | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.                                                     |
| Uniprot No.         | P24735                                                                                                            |
| Immunogen           | Recombinant Pseudomonas aeruginosa Beta-lactamase protein (27-397AA)                                              |
| Raised In           | Rabbit                                                                                                            |
| Species Reactivity  | Pseudomonas aeruginosa                                                                                            |
| Tested Applications | ELISA                                                                                                             |
| Form                | Liquid                                                                                                            |
| Conjugate           | HRP                                                                                                               |
| Storage Buffer      | Preservative: 0.03% Proclin 300<br>Constituents: 50% Glycerol, 0.01M PBS, PH 7.4                                  |
| Purification Method | >95%, Protein G purified                                                                                          |
| Isotype             | IgG                                                                                                               |
| Clonality           | Polyclonal                                                                                                        |
| Alias               | Beta-lactamase (EC 3.5.2.6) (Cephalosporinase), ampC                                                              |
| Immunogen Species   | Pseudomonas aeruginosa (strain ATCC 15692 / DSM 22644 / CIP 104116 / JCM 14847 / LMG 12228 / 1C / PRS 101 / PAO1) |
| Research Area       | Others                                                                                                            |
| Target Names        | ampC                                                                                                              |