



Caprin1 Antibody

Product Code	CSB-PA711733HA01RA
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q5M9G3
Immunogen	Recombinant Rat Caprin-1 protein (2-707AA)
Raised In	Rabbit
Species Reactivity	Rat
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Caprin-1 (Cytoplasmic activation- and proliferation-associated protein 1) (GPI-anchored protein p137) (GPI-p137) (p137GPI) (RNA granule protein 105), Caprin1, Gpiap1 Rng105
Immunogen Species	Rattus norvegicus (Rat)
Research Area	Cell Biology
Target Names	Caprin1